


# Quatro Vet Uni-Strip

***In vitro* Rapid Diagnostic Test for detection of  
Rotavirus, Cryptosporidium, Coronavirus and  
Escherichia coli K99 in bovine stool specimens**

- Specific and sensitive
- Easy to perform
- Convenient
- Simple and accurate
- No special skill for interpretation
- Minimal faecal manipulation
- No instrumentation

# Quatro Vet Uni-Strip



## Performance of the kit

### Rotavirus test:

Results on 67 faecal samples.  
Sensitivity: 92.85% - Specificity: 100%  
Reliability: 98.5%

ELISA / Rota	Positive	Negative	Total
Positive	13	0	13
Negative	1	53	54
<b>Total</b>	<b>14</b>	<b>53</b>	<b>67</b>

### Escherichia coli K99 test:

Results on 61 faecal samples  
Sensitivity : 100 % - Specificity : 100 %  
Reliability : 100 %

ELISA / E. coli K99	Positive	Negative	Total
Positive	14	0	14
Negative	0	47	47
<b>Total</b>	<b>14</b>	<b>47</b>	<b>61</b>

### Cryptosporidium test:

Results on 113 fecal samples.  
Sensitivity: 94.74% - Specificity: 92%  
Reliability: 92.9%

ELISA / Crypto	Positive	Negative	Total
Positive	36	6	42
Negative	2	69	71
<b>Total</b>	<b>38</b>	<b>75</b>	<b>113</b>

### Coronavirus test:

Results on 69 faecal samples.  
Sensitivity: 92.59% - Specificity: 100%  
Reliability: 97.10%

ELISA / Corona	Positive	Negative	Total
Positive	25	0	25
Negative	2	42	44
<b>Total</b>	<b>27</b>	<b>42</b>	<b>69</b>


For more technical information

[client.care@corisbio.com](mailto:client.care@corisbio.com)

Ordering Information

PRODUCT NAME	DESCRIPTION	STORAGE	CODE
Quatro Vet Uni-Strip	10 individual strips with 4 tests, buffers, sample taking material (test tubes, loops, tube holder)	4°C to 30°C	C-1540

— Distributed by —



CREALYS Science Park  
Rue Guillaume Fouquet 11  
5032 Gembloux - BELGIUM  
T +32 (0)81 719 910 - F +32 (0)81 719 919  
sales@corisbio.com - [www.corisbio.com](http://www.corisbio.com)



MADE IN BELGIUM

D-CUS-10-40-5/11  
March 2023